

ePURE


by  ECOLI

Fully automated nucleic acid extractor

- ❑ Liquid handling instrument based on magnetic particle extraction technology
- ❑ Purifies 1 to 12 samples simultaneously per run
- ❑ Compact bench-top design
- ❑ Built-in UV lamp decontamination
- ❑ Worry-free optimized in-built protocols
- ❑ Range of purification kits that contains all pre-filled reagents and consumables
- ❑ Covid-19 Viral RNA/DNA extraction kit
- ❑ 7" inch color touch panel
- ❑ User friendly operation, easy cleaning through front door




**READY
TO USE
FORMAT**



Cat.No.	Product name	Size	Description	Availability
E2001	ePure Blood DNA Extraction Kit	48 tests	extraction of genomic DNA from whole blood, buffy coat and leukocyte concentration	✓
E2003	ePure Viral Nucleic Acid Extraction Kit	48 tests	extraction of viral DNA/RNA from serum, plasma or cell-free fluids	✓
E2006	ePure Bacterial DNA Extraction Kit	48 tests	extraction of bacterial DNA from various sample types including environment material	on request
E2007	ePure STD Extraction Kit	48 tests	extraction of nucleic acid from cervical swab and transfer medium	✓
E2008	ePure TB DNA Extraction Kit	48 tests	extraction of MTB DNA from various sample types (e.g. sputum, BAL, CSF, urine, blood)	on request
E2004	ePure Tissue DNA Extraction Kit	48 tests	extraction of DNA from solid animal tissues, dried swab and dried blood	✓
E2009	ePure FFPE DNA Extraction Kit	48 tests	extraction of genomic DNA from formalin-fixed, paraffin-embedded tissue (FFPE) samples	on request
E2010	ePure Forensic DNA Extraction Kit	48 tests	extraction of genomic DNA from a wide range of forensic & human identity samples	✓
E2024	ePure cfDNA Extraction Kit Plus	48 tests	extraction of circulating DNA from plasma, serum or cell-free body fluids	on request

Ver.02_01/2022

 **Headquarters**
Ecoli Dx, s.r.o.
Purkyňova 74/2
11000 PRAGUE 1
Czech Republic



 **Contact**
Phone: +420 325 209 912
Cell: +420 739 802 523
E-mail: ecolidx@ecolidx.com
Web: www.ecolidx.com