



# ECOLI

## NEW LYOPHILISED KITS FOR HEPATITIS B, C AND HIV WORK SIMPLER AND FASTER!

- ◆ qualitative and quantitative detection of samples extracted from blood plasma
- ◆ internal control included
- ◆ only 2 fluorescent channels required (FAM/JOE)
- ◆ ready-to-use lyophilized reaction mixture in universal 0,2 ml tubes for quicker and easier preparation
- ◆ room temperature setup
- ◆ hepatitis kits have the same amplification profile (temperature cycling of both hepatitis kits can be performed simultaneously in one instrument)
- ◆ RUO
- ◆ 96 reactions in one kit
- ◆ variant of larger kit with 480 reactions available

### HBV

Detection limit  
3 IU/ml in a sample

### HCV

Detection limit  
7 IU/ml in a sample

### HIV

Detection limit  
25 IU/ml in a sample

#### PRODUCT NAME

#### CATALOG NUMBERS 96 REACTIONS

#### CATALOG NUMBERS 480 REACTIONS

AmpliSens® HCV-Monitor-L PCR kit. Variant FRT-L

H-4001-1-14-100-CE

H-4001-1-14-CE

AmpliSens® HBV-Monitor-L PCR kit. Variant FRT-L

H-4021-1-14-100-CE

H-4021-1-14-CE

AmpliSens® HIV-Monitor-L PCR kit. Variant FRT-L

H-4011-1-14-100-CE

H-4011-1-14-CE

**Headquarters**  
Ecoli DX, s.r.o.  
Purkyňova 74/2  
11000 PRAGUE 1  
Czech Republic



**Contact**  
Phone: +420 325 209 912  
Cell: +420 739 802 523  
E-mail: [ecolidx@ecolidx.com](mailto:ecolidx@ecolidx.com)  
Web: [www.ecolidx.com](http://www.ecolidx.com)